

Chromo-LAL

Kinetic Chromogenic Formulation

General Product Description

Chromo-LAL lysate is lyophilized with substrate reagent and buffers. It is optimized for the kinetic chromogenic LAL test method in microplate readers. Chromo-LAL is a buffered, stable and robust lysate, suitable for quantitative testing of a wide range of samples.

Sensitivity

The sensitivity for chromogenic assays is determined by the lowest standard concentration on the standard curve used for the assay. The maximum sensitivity of Chromo-LAL is 0.005 EU/mL.

Sample to Lysate Ratio

In the Chromo-LAL test, reconstituted Chromo-LAL reagent is used at a ratio of 1:1 and a volume of 100 μ L.

Performing the Test

The Chromo-LAL/sample mixture is incubated at $37 \pm 1^{\circ}$ C in a microplate reader equipped with a 405 nm filter. Software is used to construct the standard curve and calculate the endotoxin concentrations.

Reconstitution

Chromo-LAL lysate is reconstituted with LAL Reagent Water (LRW). It can also be reconstituted with Glucashield® buffer, a beta glucan inhibiting buffer, to render the assay endotoxin specific.

Stability

The reconstituted lysate is stable for 24 hours at 2–8°C or for two weeks at -20°C or colder if frozen immediately after reconstitution. Reconstituted Chromo-LAL lysate may be frozen and thawed once.

Packaging

Chromo-LAL lysate is sold by the vial for greater ordering and testing flexibility. This allows you to order the amount of lysate that is right for you.

Accessory Products

- 1. LAL Reagent Water, available in multiple package sizes, see page 4.2
- 2. Control Standard Endotoxin, *E. coli* O113:H10, 0.5 µg per vial (E0005)
- 3. Glucashield® Buffer (GB051), (1→3)-β-D-Glucan Inhibiting Buffer
- 4. Pyrotubes® (TB013), 13 x 100 mm borosilicate glass dilution tubes
- 5. Pyroplate® (CA961), 96-well microplate
- 6. 250 µL Pipette tips (PPT25), endotoxin and glucan-free
- 7. 1000 µL Pipette tips (PPT10), endotoxin and glucan-free
- 8. Incubating microplate reader

For more information please visit our website at www.acciusa.com or contact Technical Services at (800) 848–3248.

Chromo-LAL

#C0031 ----- 3.2 mL/vial

30 tests



Corporate Headquarters Associates of Cape Cod, Inc. 124 Bernard E. Saint Jean Drive East Falmouth, MA 02536-4445 T (508) 540–3444 • F (508) 540–8680 Online: www.acciusa.com UK Office
Associates of Cape Cod Int'l Inc.
Deacon Park, Moorgate Road
Knowsley, Liverpool L33 7RX, United Kingdom T (44) 151-547-7444 • F (44) 151-547-7400
Online: www.acciuk.co.uk

European Office
PYROQUANT DIAGNOSTIK GmbH
Opelstrasse 14
D-64546 Mörfelden-Walldorf, Germany
T (49) 61 05-96 10 0 • F (49) 61 05-96 10 15
Online: www.pyroguant.de

