

Product Datasheet

BVDV E2 Protein (orb1516885)

Description	Recombinant BVDV E2 Protein
Endotoxins	Not analyzed
Target	BVD Genome polyprotein
Preservatives	Lyophilized or liquid
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	Sumo-GST
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	1-140/3988
Purity	> 95% as determined by SDS-PAGE
MW	45 kDa
Uniprot ID	P19711
Expiration Date	6 months from date of receipt.