



Product Datasheet

Beta galactosidase antibody (orb100205)


Description	Beta galactosidase antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Beta galactosidase (221-290/677 aa)
Target	Beta galactosidase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	71 kDa
Uniprot ID	P16278
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Paraformaldehyde-fixed, paraffin embedde...



Paraformaldehyde-fixed, paraffin embedde...



Tissue/cell: human lung carcinoma, 4% Pa...