

Product Datasheet

BVDV E2 antibody (orb312169)

Description	BVDV E2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombined BVDV Envelope glycoprotein E2 (full length protein)
Target	BVDV E2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	P19711
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.