



Product Datasheet

Collagen-pan antibody (orb665466)

Description	Rabbit polyclonal antibody to Collagen-pan
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Collagen-pan. The exact sequence is proprietary.
Target	COL1A1;COL2A1;COL3A1;COL5A1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O88207, Q9JI03, P02461, P08121, P20908, P28481, P11087, P05539, P13941, P02458, P02454, P02452
Entrez	85490, 12842, 12831, 12825, 1280, 25412, 1281, 1277, 29393, 84032, 1289, 12824
Dilution Range	WB: 1-500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Collagen-pan ex...



Immunohistochemical analysis of Collagen...