



BOX Station One Technical Specifications

Dimensions	L 9" x W 6" x H 8.375"	
Power Supply	12V DC	
1/0	dual independent fan controls with LED indicators	
Capacity		
	Two 96 well 2D barcoded tube racks	
standard	or two microplates	
small tube accessory	52 microtubes up to 2ml capacity	
	24 mid size tubes (up to 17mm O.D.) and /or	
large tube accessory	10 large tubes (up to 30mm O.D.)	
Approx Thaw Time from -20°C (1XTE buffer)		
Full 2D barcoded tube rack		30 min
96 well microplate	7 min	
full small tube accessory	20 min	
full large tube accessory		30 min
Compatible Labware	Argos	Kimble Chase
	Axygen	Matrix
	BD Biosciences	Micronics
	Beckman Coulter	Nalgene
	Biotix	Nat. Sci. Supply
	Chemglass	Nunc
	Corning	Phenix
	Eppendorf	Sarstedt
	Falcon	Sorensen Bioscience
	Globe	Thermo Scientific
	Grenier	VWR