

***Compatible and certified with: Precision Series & Sorenson Benchtop Pipettors***



***Extends original warranty for a Lifetime!***

- Low Binding Polymer Technology reduces binding of DNA, enzymes and proteins
- Super-micron barrier inhibits pipettor and carry-over contamination
- Barrier will not trap liquid or contaminate samples if over-pipetting
- Exceptional accuracy and precision
- Sterile

**10 $\mu$ l LongReach - Working Volume: 1-20 $\mu$ l in increments of .1 $\mu$ l**

**2936T:** 10 $\mu$ l LongReach *Barrier* Pipette Tips, Low Binding, Sterile, Natural

**2945T:** 10 $\mu$ l LongReach Pipette Tips, Low Binding, Natural

**2946T:** 10 $\mu$ l LongReach Pipette Tips, Low Binding, Sterile, Natural



**250 $\mu$ l - Working Volume: 1-250 $\mu$ l in increments of 1 $\mu$ l**

**2938T:** 200 $\mu$ l NX*Barrier* Pipette Tips, Low Binding, Natural

**7244T:** 300 $\mu$ l NX Pipette Tips, Low Binding, Sterile, Natural

**7245T:** 300 $\mu$ l NX Pipette Tips, Low Binding, Sterile, Natural