

Product Datasheet

Recombinant Human HNRNPA1 Protein, N-His (orb2969321)

Catalog Number	orb2969321
Category	Proteins
Description	Recombinant Human HNRNPA1 Protein, N-His
Form/Appearance	Lyophilized
Buffer/Preservatives	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Purity	>90% as determined by SDS-PAGE.
Protein Length	Ser2-Phe372
MW	40.92 kDa
Expression System	E. coli
Endotoxins	Please contact with the lab for this information.
Biological Origin	Homo sapiens (Human)
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)