



Product Datasheet

FLOT1 antibody (orb247308)


Description	Rabbit polyclonal antibody to FLOT1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-427 of human Flotillin 1 (NP_005794.1).
Preservatives	PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	O75955
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 c...



Immunohistochemistry analysis of paraffi...



Immunohistochemistry analysis of paraffi...